



new doctoreyes ultrabright mirrors



doctoreyes presents a new generation of dental mirrors, which leaves the rest of the mirrors in the dust. With the groundbreaking ultra bright coating one can for the first time achieve close to 100% reflectivity. Even high-grade rhodium and titanium mirrors have only reached on average 75% reflectivity until now; the traditional chromium plated mirrors achieve about 65%. The mirrors come with a removable silicon handle for comfortable handling which also can be autoclaved together with the mirror.

With thorough testing the mirrors have been shaped and perfected for optimal use. Our unique mirror forms were created based on the best possible combination of reflection surface, safe handling and comfort. For making occlusal overviews we have created three different sizes.

For buccal and lingual shots we designed specially optimized shapes. For orthodontic images we offer our occlusal and buccal mirrors also in sizes which are suitable for children.



**doctorseyes
ultrabright mirrors**

doctorseyes™



occlusal mirror 12 adult size, width 65mm, glass mirror 3 mm thick, ultrabright coating,



occlusal mirror 13 adult size XL, width 70mm, glass mirror 3 mm thick, ultrabright coating,



buccal mirror 15 adult size, width 40mm, glass mirror 3 mm thick, ultrabright coating



lingual mirror 14 adult size, width 40mm, glass mirror 3 mm thick, ultrabright coating,



occlusal mirror 11 children size, width 55mm, glass mirror 3 mm thick, ultrabright coating,




occlusal mirror 02 trapezoid, 45 degree angled for mounting on doctorseye 72 and 86. mirror allows pictures up to 4 teeth. with ultrabright coating, rounded edges

occlusal mirror 04 overview, 45 degree angled for mounting on doctorseye 72 and 86, mirror allows pictures of complete upper or lower jaw.

all mirrors up to 100% reflection, rounded edges, with silicon hand grip, can be sterilized

save with our set



mirror set A consisting of buccal mirror 15, lingual mirror 14 and occlusal mirror 12, glass mirror 3 mm thick, ultrabright coating, up to 100% reflection, rounded edges, with silicon hand grip, can be sterilized



contrastor set consisting of two contrastors, different shapes, for high contrast and professional look of front teeth pictures